

THE ART AND SCIENCE OF DEFIBRILLATION™

Lifeline VIEW

The Lifeline VIEW makes it easy for anyone from first responders to untrained bystanders to respond confidently and appropriately during an emergency. This revolutionary life-saving device combines a sleek, user-friendly design with military-grade specifications—including the most proven biphasic waveform. It's also the first and only AED with a full-color interactive display that shows step-by-step videos for performing CPR, rescue breathing and external defibrillation.

You can also be confident that the device will be easy to maintain and rescue ready. The Lifeline VIEW provides clear guidance for using and maintaining the device. Additionally, it has been designed to support upgrades and enhancements.

Like the original award-winning Defibtech AED, the Lifeline VIEW includes all the mission-critical features necessary to provide the most advanced cardiac treatment and exceeds the American Heart Association performance recommendations.



martek[™]
LIFECARE
Making Life Better

Lifeline VIEW Semi-Automatic Defibrillator

Technical Specifications*

DEFIBRILLATOR

Type Semi-Automatic external defibrillator	Video Prompts Full motion video, on-screen text prompts
Model DDU-2300 Series	Controls Lighted On/Off button Lighted Shock button
Waveform Biphasic Truncated Exponential (Impedance compensated)	Display High-resolution colour LCD
Energy Adult: 150 Joules Pediatric: 50 Joules (Nominal into 50 Ohm load)	CPR Coaching Video and voice coaching On-demand video help
Charge Time* Less than 4 Seconds (Shock decision to shock)	Voice Prompts Extensive voice prompts guide user through operation of the unit
Rescue Protocol AHA 2005 Supports protocol updates by the user (password protected)	

DEFIBRILLATION/MONITORING PADS

Model Adult - DDP-2001 Child/Infant - DDP-2002	Surface Area 77 cm ² (nominal, each pad) 50 cm ² (nominal, each pad)
Type Pre-connected, single-use, non-polarized, disposable, self-adhesive electrodes with cable and connector	

BATTERY PACK

Model DBP-2003	Capacity 125 shocks or 8 hours continuous operation
Power 12V, 2800 mAh	Standby Life 4 Years (Typical, new battery at 25°C)
Type Lithium/Manganese Dioxide Disposable, recyclable, non-rechargeable	Low Battery Indicators Visible, Audible

PHYSICAL

Size 7.3 x 9.5 x 2.3 inches (18.5 x 24 x 5.8 cm)	Weight (Approximate) Less than 3 lbs (1.4kg) (With battery)
---	--

ENVIRONMENTAL

Relative/Humidity Operating/Standby: 5%- 95% (non-considering)	Altitude -500 to 15,000 ft (-150 to 4500m) per MIL-STD-810F 500.4 Procedure II	Sealing/Water Resistance IEC60529 class IP54; Splash proof, Dust & waterjets Protected (Battery Pack installed)	EMC (Immunity) IEC 61000-4-3 and IEC 61000-4-8
Temperature Operating: 0 to 50°C (32 to 122°F)	ESD EN 61000-4-2(open air up to 8kV or direct contact up to 6kV)	Crush Test 1,000 pounds (450Kg)	EMC (Emission) CISPR 11 Group 1 Level B and FCC Part 15
Shock/Drop Abuse Tolerance MIL-STD-810F 516.5 Procedure IV (1.2 meter, any edge, corner or surface, in standby mode)	Vibration Ground (MIL-STD-810F 514.4 Category 20) Helicopter (RTCA/DO-160D, Section 8.8.2, Cat R, Zone, 2, Curve G) Jet Aircraft (RTCA/DO-160D Section 8. Cat H, Zone 2, Curves B & R)		

SELF TESTS

Automatic Automatic daily, weekly and monthly circuitry tests
Battery Insertion System integrity test on battery insertion
Pad Presence Pads preconnected tested daily
User-Initiated Unit and battery pack system test may also be initiated by the user
Status Indication Visual and audible indication of unit status
Status Screen Unit self-test results Pad and battery information (status and expiration)

EVENT DOCUMENTATION

Internal Event Record Critical ECG segments and rescue event parameters are recorded and can be downloaded to a removable data card
PC-Based Event Review ECG with event tag display, and audio playback when available
Removable Storage (Optional) Up to 30 hours of ECG and event data storage (no audio option) or up to 3 hours of audio, ECG and event storage (audio option) on a removable data card. Actual length of storage is dependant on card capacity.
USB Port Event download and maintenance operations

PATIENT ANALYSIS SYSTEM

Patient Analysis Automatically evaluates patient impedance for proper pad contact. Monitors signal quality and analyzes patient ECG for shockable/non-shockable rhythms
Sensitivity/Specificity Meets AAMI-DF-80 specifications and AHA recommendations